

AudBase – ECAA Compatible Equipment

Instrument	Description	Connection/Data Transfer Type	PC Controlled
GSI 61	Audiometer	**Serial (USB Converter)	No
GSI AudioStarPro	Audiometer	XML	No
GSI Pello	Audiometer	XML	No
GSI Tymptstar	Tympanometer	**Serial (USB Converter)	No
GSI TymptstarPro	Tympanometer	XML	No
GSI 1739	Aud/Tymp Combination	USB	No
Interacoustics AC40 (old style)	Audiometer	**Serial (USB Converter)	No
Interacoustics AC40 (new style)	Audiometer	XML	No
Interacoustics Equinox	Audiometer	XML	*Yes
Interacoustics Affinity	Audiometer	XML	*Yes
Interacoustics Titan	Tympanometer OAE	XML	No
Maico Easy Tymp	Tympanometer	XML	No
Otometrics Astera	Audiometer	XML	*Yes
Otometrics OtoFlex	Tympanometer	XML	No
Otometrics Zodiac	Tympanometer	XML	*Yes
Otometrics Capella	OAE	XML	*Yes
Otometrics Aurical (via Otsuite)	Audiometer	XML	*Yes
Otometrics Conera (via Otsuite)	Audiometer	XML	*Yes

*Instruments that are PC controlled will require the device application to be loaded on PC.

Revised 2020

**Serial instruments will use a serial to USB converter and have COM ports configured.

Contact us:

1-888-959-6790

support@audbase.com

www.audbase.com

